Biomass Pyrolysis in a Heated Microtubular Reactor

AnGayle K. Vasiliou,^{1,3} Krzysztof M. Piech,¹ Beth Reed,¹ Xu Zhang,² Mark R. Nimlos,³ Musahid Ahmed,⁴ Amir Golan,⁴ Oleg Kostko,⁴ David L. Osborn,⁵ John W. Daily,⁶ John F. Stanton,⁷ and <u>G. Barney</u> Ellison* ¹

- Department of Chemistry & Biochemistry University of Colorado Boulder, CO 80309-0215
- * Corresponding author: <u>barney@jila.colorado.edu</u>
- Jet Propulsion Laboratory California Institute of Technology 4800 Oak Grove Drive Pasadena, CA 91109-8099
- National Renewable Energy Laboratory 1617 Cole Blvd.
 Golden, CO 80401
- Chemical Sciences Division LBNL MS 6R-2100 Berkeley CA 94720

- Combustion Research Facility Sandia National Laboratories PO Box 969 MS 9055 Livermore, CA 94551-0969
- ⁶ Center for Combustion & Env. Research Department of Mechanical Engineering University of Colorado Boulder, CO 80309-0427
- Institute for Theoretical Chemistry Dept. Chemistry, University of Texas Austin, TX 78712

A heated SiC microtubular reactor has been developed to decompose biomass monomers such as CH₃CHO, C₆H₅OH, C₆H₅OCH₃, HOC₆H₄OCH₃, C₆H₅CH₂CH₂OC₆H₅, and furan. The pyrolysis experiments are carried out by passing a dilute mixture of the organic substrate (roughly 0.1-1%) entrained in a stream of a buffer gas (either He or Ar) through a heated SiC reactor that is 2-3 cm long and 1 mm in diam. Common pressures in the reactor are 50-200 Torr and the SiC tube is heated in the range of 1200-1900 K. Typical transit times through the reactor are 50-200 μ sec after which the gas mixture emerges as a skimmed molecular beam at a pressure of approximately 10 μ Torr and all chemistry is quenched. The reactor has been deployed in pulsed and CW configurations. In Colorado a pulsed reactor is used because the pyrolysis products are

identified by photoionization mass spectroscopy with a 10 Hz YAG laser at $\lambda_0 = 118.2$ nm or 10.487 eV. Separate experiments use matrix absorption infrared spectroscopy to identify the pyrolysis products and to confirm the assignments of the PIMS. In Calif. a CW reactor is used because a CW synchrotron is used as the light source for the PIMS. The pyrolysis of CH₃CHO will be discussed. observe CH₃CHO (+M) → $CH_3 + H + CO + CH_2 = C = O$, CH_2 =CHOH, HC= $CH + H_2O$.

