Thermodenuder-AMS, and Organic Component Analysis, and CCN Closure with AMS Data

 Peter DeCarlo, Alex Huffman, Qi Zhang, James Allan, Mike Cubison, Eben Cross, Doug Worsnop, Paul Ziemann, Athanasios Nenes, Anne Jefferson, Betsy Andrews, John Ogren, and Jose-Luis Jimenez

Univ. of Colorado-Boulder, Univ. of Manchester, Aerodyne Research, UC-Riverside, Georgia Tech, NOAA CMDL

Chebogue Pt. Data Analysis Mtg – NOAA 3/10/2005

Part I: Thermodenuder-AMS

- Idea: measure change in aerosol composition as a function of volatilization temperature
- Focus on organics, "polymers"?

Themodenuder



Design after Wehner et al. J. Aerosol Sci. 33: 1087-1093.(2002)

Thermodenuder Temperature Profile



Thermodenuder/TDPBMS Monocarboxylic and Dicarboxylic Acids



TD-AMS Time Series



Appears consistent with VTDMA and TAG/Org fraction

Part II: Organic Component Analysis

- Some chemical resolution on the whole organic aerosol
- •Apply method developed for urban areas
- "Custom" Principal Component Analysis
- Qi Zhang applying this to 25+ worldwide datasets

OOA HOA Quantification Algorithm (see Poster 12PA-11)



HOA and OOA Time Series in <u>Pittsburgh</u>



Chebogue: Results of First Pass



SO₄/Org vs Org Components



Part III: Aerosol-CCN Closure

- Apply method developed for urban areas
- "Custom" Principal Component Analysis
- Qi Zhang applying this to 25+ worldwide datasets

Data Integration Algorithm







Closure Strategy







Period 2

Lots of organics.

mode 2

Number = 1374.1 Diameter = 70 Sigma = 1.64MF Solubles: 0.13413 MF Slightly Soluble: 0.82963 mode 3 Number = 237.17 Diameter = 297.97 Sigma = 1.5342 MF Solubles: 0.10462 MF Slightly Soluble: 0.85877



Mode 1	Mode 2	Mode 3
Col. Eff. 1.00 🔶	0.35 🔶	0.34 🔶
EC / Pri 0.10 🕀	0.30 🕀	0.00
Density [lon. Eff.]	1.15 () 1.40 ()	1.20 ↔ 1.40 ↔

lon. Eff. 1.40 🕀

SuICE 1.00 🕀
Nit CE 1.00
AmmCE 1.00
OrgCE 1.00





Period 5 Lots of Sulfate

mode 2			
Number 0745 4 (em2)			
Number = 2715.1 /cm3			
Diameter = 105.34 nm			
Sigma = 1.4174			
ME Solubles: 0.58821 ME Slightly Soluble: 0.40131			
mode 2			
Number = 1132.6 /cm3			
Diameter = 207.13 nm			
Sigma = 1.5033			
MF Solubles: 0.65966 MF Slightly Soluble: 0.33014			
Lens Transmission Curve Liu 💌			
Weighing Type for Fit Standard 💌			
☑ Weigh for Pt Spacing			
User Defined Weighing Factors			
(Use when weighing is set to "Defined")			
Number 1.0 🔂 Volume 1.0 🕀			
Sulfate 1.0 🕀 Nitrate 1.0			
Ammonium 1.0 🕀 Organic 1.0 🕀			
Miscellaneous Parameters			

Mode 1	Mode 2	Mode 3
Col. Eff. 1.00 ↔	0.57 🕀	0.63 ↔
EC / Pri 0.10 ↔	0.30 🕀	0.00 ↔
Density [1.15 () 1.40 (⊖	1.20

SuiCE 1.00 🕀 Nit CE 1.00 🕀 AmmCE 1.00 🕀 OrgCE 1.00 🕀



Conclusions + Future Work

- TD-AMS
 - So far appears consistent with trends of study
 - Next: full analysis of magnitude & chemistry of volatiliz.
 - Compare to VTDMA and TAG/Org fraction
- Organic component analysis
 - Mostly OOA, two types, correlated with other obs.
 - Developing procedures to separate 3rd component
- CCN closure
 - Reasonable results, but lots of "knobs"
 - Next: Mike Cubison, OPC size dist
- Light scattering closure
- Hygroscopicity closure